



MicroSnap™

Coliform



MicroSnap Coliform is a rapid test for detection and enumeration of Coliform bacteria. The test uses a novel bioluminogenic reaction that generates light when enzymes that are characteristic of Coliform bacteria react with specialized substrates. The light generating signal is then quantified in the EnSURE luminometer. Results are available in 8 hours or less, depending upon required level of detection. Single figure organisms can be detected in 8 hours, delivering results in the same working day or shift.

Benefits:


- 6 hour quantitative results (CFU counts)
- 8 hour presence/absence results
- Broad CFU detection range (10 – 10,000 CFU) means fewer dilutions and less sample prep saving time and money
- Suitable for product or environmental samples
- Pass/Fail result at required detection levels are easy to understand
- Uses common microbiology sample preparation for product samples
- Same day results allow for timely corrective action
- In-house testing reduces costs from outside lab fees
- AOAC Research Institute *Performance Tested Method* (#071302)
- 12 month shelf life at refrigerated temperatures (2-8°C)

Sensitivity & Specificity:

- Detects down to <10 CFU in 8 hours
- Detects Coliform bacteria including *Escherichia*, *Klebsiella*, *Citrobacter*, and *Enterobacter*

CFU : RLU - Conversion

CFU/ml or g		Coliform (RLU)
<10	≈	<2
<20	≈	<4
<50	≈	<7
<100	≈	<12
<200	≈	<20
<500	≈	<35
<1,000	≈	<60
<5,000	≈	<180
<10,000	≈	<300



YouTube [youtube.com/HygienaTV](https://www.youtube.com/HygienaTV)

Watch an instructional demo at [youtube.com/HygienaTV](https://www.youtube.com/HygienaTV)



Find support documents, instructional videos, and more at www.hygiena.com


keydiagnostics

T: 02 8212 4074 F: 02 9423 6992
info@keydiagnostics.com.au
www.keydiagnostics.com.au
 PO Box 1038, Gympie, NSW, 2227

Catalog No.	Description	Quantity
MS-E-CEC	MicroSnap Coliform & <i>E. coli</i> Enrichment Device	100
MS-E-CEC-BROTH	MicroSnap Coliform & <i>E. coli</i> Enrichment Broth	100
MS-COLIFORM	MicroSnap Coliform Detection Device	100

REV B 072014

MicroSnap Coliform Procedure: Watch an instructional video at [youtube.com/HygienaTV](https://www.youtube.com/HygienaTV)



Related Products:



EnSURE

EnSURE is a high-sensitivity luminometer that measures multiple quality tests such as ATP, TVC, Coliform, *E. coli*, and more all on one system. The system includes SureTrend data tracking software to record, trend and analyze test results.

MicroSnap *E. coli*

MicroSnap *E. coli* is a rapid test for detection and enumeration of *E. coli* bacteria. Results are available in 8 hours or less. MicroSnap *E. coli* detects down to 1-10 CFU of *Escherichia* bacteria.

MicroSnap Total

MicroSnap Total delivers total viable organisms test results in 7 hours, giving you the advantage of knowing microbial levels in the same working day or shift. MicroSnap Total detects down to <100 CFU of Gram positive and Gram negative aerobic and facultative anaerobic bacteria.

MicroSnap Enterobacteriaceae (EB)

MicroSnap EB delivers TVC test results in 7 hours, giving you the advantage of knowing microbial levels in the same working day or shift. MicroSnap EB detects down to 50 CFU Enterobacteriaceae bacteria and can be used as an indicator for *Salmonella* spp.

Catalog No.	Description	Quantity
ENSURE	EnSURE Monitoring System with SureTrend Software	1
See catalog	MicroSnap <i>E. coli</i>	100
See catalog	MicroSnap Total	100
See catalog	MicroSnap Enterobacteriaceae	100