

Get truly representative liquid samples as easily as 1,2,3

Sanitary Ports
QualiTru stainless steel sanitary ports are manufactured under rigorous quality control standards and allow easy access for aseptic process monitoring.

1



2

Sterile Septa
Our sterile septa allow aseptic samples to be drawn a single time from each channel virtually eliminating the risk of contamination.

Sterile Collection Units
We offer a complete line of single-use, sterile disposable collection units in a variety of sizes and configurations for every application.

3



ASEPTIC SAMPLING SYSTEM
Complete Product Guide



651-501-2337
sales@qualitru.com
qualitru.com



3-A Sanitary Standard Certified

QualiTru Sampling Systems regulatory approvals include:
FDA M-I-06-12 Needle Only Method • FDA M-I-16-17 Milk Tank Trucks • FDA M-I-06-6 Direct Load Tankers

keydiagnostics

T: 02 8212 4074 F: 02 9423 6992
info@keydiagnostics.com.au
www.keydiagnostics.com.au
PO Box 1038, Gymea, NSW, 2227

TRUSTREAM

TruStream7 Sanitary Ports

TruStream7 Tank Port	212113
TruStream7 Tri-Clamp Port 1.5"	211125
TruStream7 Tri-Clamp Port 2"	211129
TruStream7 Tri-Clamp Port 2.5"	211133
TruStream7 Tri-Clamp Port 3"	211137
TruStream7 Tri-Clamp Port 4"	211145

TruStream7 Tri-Clamp Elbow 1"	213019
TruStream7 Tri-Clamp Elbow 1.5"	213025
TruStream7 Tri-Clamp Elbow 2"	213029
TruStream7 Tri-Clamp Elbow 2.5"	213033
TruStream7 Tri-Clamp Elbow 3"	213037
TruStream7 Tri-Clamp Elbow 4"	213045

TruStream7 Recessed Port 4" x 2"	212123
TruStream7 Recessed Port 4" x 4"	212127
TruStream7 Recessed Port 5.5" x 4"	212133
TruStream7 Recessed Port 5.5" x 5"	212135

TruStream7 Recessed 5.5" x 5" w/Security Lid	212433
--	--------

TruStream7 Butt Weld Elbow 1.5"	213075
TruStream7 Butt Weld Elbow 2"	213079
TruStream7 Butt Weld Elbow 2.5"	213083
TruStream7 Butt Weld Elbow 3"	213087
TruStream7 Butt Weld Elbow 4"	213095

TruStream7 Tri-Clamp Tee 1"	215141
TruStream7 Tri-Clamp Tee 1.5"	215145
TruStream7 Tri-Clamp Tee 2"	215147
TruStream7 Tri-Clamp Tee 3"	215151
TruStream7 Tri-Clamp Short Tee 1.5"	215515
TruStream7 Tri-Clamp Short Tee 2"	215520

TruStream7 Butt Weld Tee 1"	215071
TruStream7 Butt Weld Tee 1.5"	215075
TruStream7 Butt Weld Tee 2"	215077
TruStream7 Butt Weld Tee 2.5"	215079
TruStream7 Butt Weld Tee 3"	215085
TruStream7 Butt Weld Tee 4"	215089

TruStream7 Wrench	310011
TruStream7 Stainless Steel Nut	219021
TruStream7 Dust Cover	310121

TruStream12 Sanitary Ports

TruStream12 Tank Port	212214
TruStream12 Tri-Clamp Port 2"	211164
TruStream12 Tri-Clamp Port 2.5"	211168
TruStream12 Tri-Clamp Port 3"	211174
TruStream12 Tri-Clamp Port 4"	211184

TruStream12 Recessed Port 6" x 2"	212224
TruStream12 Recessed Port 6" x 3"	212228

TruStream12 Tri-Clamp Elbow 2"	213124
TruStream12 Tri-Clamp Elbow 2.5"	213128
TruStream12 Tri-Clamp Elbow 3"	213134
TruStream12 Butt Weld Elbow 2"	213164
TruStream12 Butt Weld Elbow 2.5"	213168
TruStream12 Butt Weld Elbow 3"	213174

TruStream12 Wrench	310022
TruStream12 Stainless Steel Nut	219054
TruStream12 Dust Cover	310162

TruStream Sanitary Ports

Our 3-A Certified sanitary ports are designed to work in a wide variety of locations including processing lines, silos, tanks and milk tanker trucks. QualiTru sanitary ports can also be customized to your specification.

TruStream7
Tri-Clamp Elbow 2" #213029

TruStream12
Tri-Clamp Port 2" #211164



TruStream12 Septa #110084

TRUSTREAM12



Place your order today
651-501-2337
sales@qualitru.com
qualitru.com



TruStream7 Septa #110011

TRUSTREAM7

TruStream
12 Channel Sterile Septa

Our 12-channel sterile septa are built for durability and ease-of-use. Each single-piece septa is made of food-grade materials and features a permanent, thermal-aluminum label that clearly indicates which channels have been used. Once all 12 channels have been used, simply replace with a new sterile septa. There is no better way to gather aseptic samples.

TruStream
7 Channel Sterile Septa

Simple to install. Simple to use. Simple to see when a channel has been used. The TruStream7 channel sterile septa is made from food-grade contact materials and allows aseptic samples to be drawn a single time from each channel. Once all 7 channels have been used, simply replace with a new sterile septa.

TruStream7 Sterile Septa

TruStream7 Septa	110011
TruStream7 Septa (high temp)	110021

New TruStream12 Sterile Septa

TruStream12 Sterile Septa

TruStream12 Septa	110064
TruStream12 Septa (high temp)	110074



TruStream7 Syringe Kit #114108

TRUSTREAM7

QualiTru Sterile Collection Units

TruDraw	Single Sampler 2oz	112021
---------	--------------------------	--------

TruStream	Sampling Assembly 250ml	111450
TruStream	Sampling Assembly 2L	111730
TruStream	Sampling Assembly 7L	111930

TruMotion	Sampling Assembly 2L w/2.0mm silicone tubing for Peristaltic Pump use.....	111750
TruMotion	Sampling Assembly 2L w/3.2mm silicone tubing for Peristaltic Pump use.....	111770

TruStream7	Syringe Kit.....	114108
------------	------------------	--------

Other Consumables

QualiTru	Needles 14g	113314
QualiTru	Needles 16g	113316
QualiTru	Needles 18g	113318
QualiTru	Needles 20g	113320
QualiTru	Needles 22g	113322
QualiTru	Needles 23g	113323

QualiTru	Tamper Evident Device	310021
QualiTru	Security Seals	310020
QualiTru	Flow Clamp	112010
QualiTru	Flow Clamp Valve.....	112035

QualiTru Sterile
Collection Units

QualiTru offers a complete line of sterile, disposable sample bags in a variety of sizes making them ideal for all your sampling needs. All our collection units are built tough and provide a worry free handling, storing and shipping solution.

TruDraw Single Sampler 2oz. #112021

TRUDRAW
SINGLE SAMPLER

The TruDraw Single Sampler is a sterile, sealed container with an attached needle. Our new TruDraw Single Sampler is leak-proof and includes security features and ease-of-use for drawing samples up to two ounces. This industry first sampling method provides sample tracking capabilities and virtually eliminates the risk of external contamination. For added security, each TruDraw has a tamper evident feature. These simple to use samplers work with all QualiTru TruStream Septa.



TruStream
Collection Bags

TruStream sterile collection bags come in a variety of sizes and configurations for any sampling need. Each bag is equipped with a needle and tubing and work with all QualiTru TruStream Septa.

TRUSTREAM

TruStream Sampling Assembly 2L
#111730

