

Siltex Pratika S20: Laboratory Sterilizer

New to Siltex, Pratika's qualifications are impressive: 20 litre chamber, Type S Sterilizer that is suitable for liquids, SD card for electronic cycle archiving, clean water for each cycle, an external waste water tank, optional printer and user determined cycle times.



- Type S sterilizing
- Capacity for 10 x 250ml Schott Bottles
- Capacity for 4 x 500ml Schott Bottles
- Capacity for 2 x 1litre Schott Bottles (Optional Tray)
- 121°C or 134°C sterilizing cycles
- Closed door drying cycle
- Digital temperature readout
- SD card writer and SD card
- User selectable cycle lengths
- Printer option available
- External dump tank
- 20 Litre Chamber
- Optional 60minute controller for certain applications where liquid or agar is used

Sterilizer comes with circular rack system and 4 trays, tray handle, 2 wire racks, SD card and SD card reader, dump tank and technical manual. Printers require one additional power point, and connect directly to the sterilizer.

Technical specifications:

Unit dimensions: 45 x 39 x 58cm (w x h x d) Unit weight: 45kg
 Chamber capacity: 20 Litre - 24cm diameter, 43 cm length
 5 litre external dump tank Standard 240V, 10A power outlet required



T: 02 8212 4074 F: 02 9423 6992
www.keydiagnostics.com.au
 E: infoxkeydiagnostics.com.au
 PO Box 1038 Gympie NSW 2227