

Water activity accessories

Disposable Sample Cups

Capacity: 15mL; Fill: 7mL or less recommended.
Material: High Density Polyethylene (HDPE).
Lids: allow the short term storage of samples.

You receive 100 sample cups with your order. Additional sample cups can be purchased in boxes of 500 with discounts for multiple box orders.



The sample cups are intended to be single-use disposable cups. We cannot recommend that you wash and reuse your sample cups. There is a chance a washed cup may still have residual contamination from the previous sample or may not be completely dry from washing.

Water Activity Standards

When you own AquaLab, you own the fastest water activity meter in the industry. To ensure that your AquaLab stays that way, you can check its performance with our ready-to-use Performance Verification Standards. You no longer need to purchase and store reagent grade salts, or own a high-precision balance.

Verification Standards are available in four water activity levels: 0.984, 0.760, 0.500, and 0.250 aw.

Avoid handling and mixing errors.

AquaLab customers don't like mixing their own salt standards. They complain it's messy and time-consuming. They also question the quality of their solutions. Now, you can verify AquaLab's performance with confidence. AquaLab Performance Verification Standards are produced under a strict quality assurance regime. An independent third party verifies the accuracy of each solution.

Precision salt solutions offer low-cost quality insurance for you.

Using AquaLab is easier than ever. Prepackaged salt standard solutions are immediately accessible for performance verification, saving you time and money. Operate your instrument with certainty and insure the quality of your product.



■ FOR FURTHER INFORMATION CONTACT:



T: 02 8212 4074 F: 02 9423 6992
info@keydiagnostics.com.au
www.keydiagnostics.com.au
PO Box 1038, Gympie, NSW, 2227